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English Language

1. Discription

The high purity manual valve is made from a stainless steel bonnet and a stainless steel hand-wheel. The design is following the GMP demands for sterile applications in consequent way.

Excellent use for sterilisation processes with a fingerless compressor with a defined sealing line between diaphragm, actuator bonnet and valve body. It is a design which provides high efficiency for sterilisation processes and a reliable sealing with minimum bolt torque.

Smooth outside surfaces and no explored nuts and other screw parts guaranty an optimised wash down design in every assembly position.

The sealed manual actuator has an integrated optical indicator not disturbing the operation of the hand-wheel. The valve bonnet is sealed before the spindle and the hand-wheel is a non rising.

1.1 Technical Data:

| | |
|----------------------|---|
| Actuation manual: | Hand operated actuator with stainless steel bonnet and stainless steel hand-wheel non rising By turning the hand-wheel clockwise the valve is closing By turning the hand-wheel anti- clockwise the hand- wheel is opening. |
| Working pressure: | 0 – 10 bar |
| Working temperature: | 160°C depending on diaphragm |
| Diaphragm material: | EPDM; FPM; PTFE |
| Valve body: | Stainless Steel Investment Cast 1.4435 (316L) Stainless Steel forged 1.4435 (316L) From solid block (various materials) |
| End connection: | see TD 970009 |

Directions

The legal rules and guidelines have to be followed by choosing the proper materials

1. General

All valves are tested before leaving the factory to ensure full specified performance and bubble tight shut off. The valves are packed according internal specification if there are no special requests from the customer. Before installing the valves make sure visually there are no damages.

2.1 Installation

The valve can be installed in any position and direction. If drain-ability is an issue there are data sheets or drawings on request available.

Before installing make sure the applied materials for the working media and considering the temperature are suitable for the application.

Depending on the overall Conditions and the weight of the valve a reasonable fastening has to be done.

All valves are marked according EN 19.

For technical description, dimensions and the min. and max. working pressure with temperature see data sheet TD 970009 and TD 060004.

After start up the plant, make sure the diaphragm valve does not leak between body and diaphragm eventually re fix the bolt connection (12).

4. Diaphragm change

4.1 Security Directions

If the valve is already installed or in line, intercept the conveyed fluid upstream of the valve and ensure it is not under pressure. If necessary relax the system and drain downstream in the proper place. If temperature is applied take care the valve and the system is cooled down under the evaporation temperature of the media to avoid scalds. In addition at poisonous or aggressive media the valve has to be ventilated.

4.2 Disassembly

The diaphragm valve is the most forced part in the diaphragm valve. The stress and wear is mechanically or chemically caused by the media. The cycles of checking the diaphragm therefore should be ruled depending on the working conditions and the cycle time of the valve itself. The check of the diaphragm can be done by disassembling the actuation from the body.

Unscrew the four bolts (12) in order to separate the body (1) from the actuator and unscrew the diaphragm (2) from the stem. If needed clean or exchange the diaphragm (2) see the assembly description.

4.3 Assembly

The diaphragm (2) should be screwed in clockwise until resistance is felt upon which the diaphragm should be screwed anti-clockwise until alignment of the bolt hole centre (12) is achieved.

Fix the manual actuator with the screws (12) onto the body.

Tighten the bolts (12) cross over wise and make sure the diaphragm is not over pressed.

Tighten the bolts with a torque key according the recommended torques showed in table (1)

| Tabelle 1/ Table 1 Nennweite Size (mm) | Torques for the screws (12) in Nm Anzugsdrehmomente für Ventilkörperschrauben (12) in Nm | |
|--|---|--|
| | Membranwerkstoff Diaphragm material PTFE or TFM | Membranwerkstoff Diaphragm material Elastomere e.g. EPDM |
| 397.15 | 2 | 1,5 |
| 997.15 – 25 | 6 | 5 |
| 997.32 - 40 | 12 | 9 |
| 997.50 | 18 | 14 |

5. Accessories


As option an adjustable travel stop or/ and an adjustable stroke limiter is available. The parts are integrated in the hand-wheel and not visible which guaranties all the advantages remain as far as the design for high purity applications are concerned.

6. Maintenance

SED diaphragm valves are maintenance free.

For all work on the valve the security directions Point 4.1 has to be strictly followed.

If other languages for the above instructions are requested advise.

| Typ | A | $\varnothing B$ | C | L |
|--------|-----|-----------------|----|-----|
| 997.15 | 103 | 92 | 19 | 120 |
| 997.20 | 103 | 92 | 19 | 120 |
| 997.25 | 103 | 92 | 19 | 120 |
| 997.32 | 140 | 135 | 26 | 153 |
| 997.40 | 140 | 135 | 26 | 153 |
| 997.50 | 150 | 135 | 32 | 173 |
| 397.15 | 77 | 62 | 12 | 108 |